

FELDSPATHIC VENEERS: THE GOLD STANDARD IN DENTAL AESTHETICS

Manav Kalra

In an era dominated by digital workflows and high-strength ceramics, feldspathic porcelain remains the benchmark against which all aesthetic restorative materials are measured. While materials such as lithium disilicate (e.max) have gained widespread popularity due to their strength and versatility, when the primary objective is the creation of the most natural and life-like smile possible, feldspathic porcelain continues to be the gold standard.

The unique advantage of feldspathic ceramic lies in its exceptional optical properties. Its ability to mimic the translucency, fluorescence, opalescence, and depth of natural enamel is unmatched by any pressed or milled ceramic. Because feldspathic veneers are meticulously layered by skilled ceramists, every restoration can be individually characterised to replicate the subtle nuances found in natural teeth.

Contrary to common perception, strength is rarely the limiting factor in veneer dentistry. Modern adhesive protocols have transformed thin feldspathic restorations into highly predictable and durable restorations when bonded to enamel. Numerous long-term clinical studies have demonstrated excellent survival rates, particularly when case selection and adhesive techniques are respected.

Lithium disilicate veneers undoubtedly offer greater flexural strength and are invaluable in situations requiring additional masking or structural support. However, this strength comes at the expense of some optical refinement. For the discerning clinician pursuing the highest level of aesthetics, feldspathic porcelain provides a level of vitality and natural beauty that remains difficult to surpass.

What is particularly exciting today is the remarkable level of craftsmanship being demonstrated by dental ceramists across India. A growing number of highly skilled technicians are embracing traditional feldspathic layering techniques, working on refractory dies and utilising Geller models to achieve exceptional precision and artistry. These analogue techniques demand patience, experience, and a deep understanding of dental morphology, yet they allow ceramists to create restorations with a level of individualisation and natural beauty that is difficult to replicate through purely digital workflows. Their dedication is helping elevate the standard of aesthetic dentistry in India to a level that rivals the finest work being produced anywhere in the world.

The following case demonstrates the power of minimally invasive dentistry combined with masterfully crafted feldspathic porcelain veneers—showing why, decades after its introduction, this material continues to define excellence in aesthetic restorative dentistry. ■

Lab Credits: Katara Dental Pvt Ltd

Ceramist: Mr. Vasin



Fig 1: Pre-operative extraoral view.



Fig 2: Pre-operative intraoral view showing existing restorations and aesthetic concerns.




Fig 3: Cross-polarization image revealing the underlying restoration and substrate characteristics.



Fig 4: Cross-polarization photography used for accurate shade analysis and ceramic communication.



Fig 5: Rubber dam isolation during tooth preparation enhances visibility and soft tissue control.




Fig 6: Rubber dam retraction facilitates precise subgingival margin preparation.



Fig 7: Conservative veneer preparation design.



Fig 8: Interproximal design: Mesial slice & distal wing preparation



Fig 9: More preparation views



Fig 10: Retraction cord placed to achieve optimal soft tissue displacement before the impression.



Fig 11: Detailed view of the gingival retraction prior to impression making.



Fig 12: Hand-layered feldspathic veneers fabricated on Geller models.




Fig 13: Feldspathic porcelain veneers demonstrating precise morphology and surface texture.



Fig 14: Completed feldspathic restorations on Geller models prior to delivery.



Fig 15: Intraoral veneer try-in using water-soluble trial paste.



Fig 16: Rubber dam isolation and air-abrasion protocol prior to adhesive bonding.

ARTICLE CITATION

XX



Fig 17: Immediate post-operative result following veneer cementation.



Fig 18: Two-week post-operative review showing soft tissue maturation.



Fig 19: Occlusal view demonstrating integration from palatal side



Fig 20: Two-week post-operative smile view.



Fig 21: Pre- and post-operative intraoral comparison.



Fig 22: Pre- and post-operative extraoral comparison.



Fig 23: Final smile showcasing the aesthetic integration of feldspathic porcelain veneers.

ABOUT THE AUTHOR



Dr. Manav Kalra is a private practitioner and dental educator based in New Delhi. He completed his postgraduate training in Prosthodontics and Oral Implantology from Manipal University. His clinical interests include Aesthetic Dentistry, Full Mouth Rehabilitation, and Implantology, with a focus on minimally invasive, additive techniques.

Dr. Kalra is the founder of the Kalra Dental Academy, where he mentors participants through his signature course, The Art and Science of Porcelain Veneers. He is also associated with 3M as a speaker and regularly shares his knowledge through lectures and workshops at institutions across India.